

Product Datasheet

IFN-gamma receptor 1 antibody (orb1089965)

Description	Rat monoclonal antibody to IFN-gamma receptor 1
Species/Host	Rat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	Blocking, FC
Immunogen	This antibody was raised by immunising rats with interferon gamma receptor.
Target	IFN-gamma receptor 1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	This antibody is recommended as an efficient blocking agent for the interferon gamma receptor (Basu et al. 1988). It has been also used in the flow cytometric analyses of this receptor protein.
Isotype	IgG2a
Clonality	Monoclonal
Purity	Purified
Uniprot ID	P15261
Expiration Date	12 months from date of receipt.