

## Product Datasheet

### IFN-gamma receptor 1 antibody (orb1089966)

<b>Description</b>	Rabbit monoclonal antibody to IFN-gamma receptor 1
<b>Species/Host</b>	Rat
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	Blocking, FC
<b>Immunogen</b>	This antibody was raised by immunising rats with interferon gamma receptor.
<b>Target</b>	IFN-gamma receptor 1
<b>Preservatives</b>	PBS with 0.02% Proclin 300.
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
<b>Note</b>	For research use only
<b>Application notes</b>	This antibody is recommended as an efficient blocking agent for the interferon gamma receptor (Basu et al. 1988). It has been also used in the flow cytometric analyses of this receptor protein.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified
<b>Uniprot ID</b>	<a href="#">P15261</a>
<b>Expiration Date</b>	12 months from date of receipt.