

Product Datasheet

SNCG antibody (orb1089985)

Description	Rabbit monoclonal antibody to SNCG
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	The original antibody was raised against the KLH-conjugated peptide "TVEEAENIAVTSGVVR", derived from Gamma-synuclein, which plays a role in neurofilament network integrity.
Target	SNCG
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	The original characterization data for this antibody can be found here .
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q6FHG5
Expiration Date	12 months from date of receipt.