

Product Datasheet

XPNPEP1 antibody (orb1089998)

Description	Rabbit monoclonal antibody to XPNPEP1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Peptide "TSLSLDEVYLIDSGAQYK" derived from Xaa-Pro aminopeptidase 1 conjugated to KLH.
Target	XPNPEP1
Preservatives	PBS with 0.02% Proclin 300.
Concentration	1 mg/ml
Storage	Store at 4°C for up to 3 months. For longer storage, aliquot and store at -20°C.
Note	For research use only
Application notes	Original data characterizing this antibody may be found here .
Isotype	IgG
Clonality	Monoclonal
Purity	Purified
Uniprot ID	Q9NQW7
Expiration Date	12 months from date of receipt.