


Product Datasheet

Notum antibody (orb109018)

Description	Rabbit polyclonal to Notum
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Preservatives	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Form/Appearance	Provided as solution in phosphate buffered saline with 0.08% sodium azide
Concentration	See vial for concentration
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	Antibody can be used for Western blotting (5-10 µg/ml) and immunohistochemistry. Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG
Source	Rabbits were immunized with recombinant full-length notum protein from mouse.
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cell lysa...