



Product Datasheet

Skp2 antibody (orb10902)


Description	Rabbit polyclonal antibody to Skp2
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, WB
Immunogen	Recombinant fusion protein of human SKP2 (NP_005974.2).
Preservatives	PBS with 0.05% proclin300,50% glycerol,pH7.3.
Storage	Store at -20°C. Avoid freeze / thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48KDa
Uniprot ID	Q13309
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U2OS cell...



Immunohistochemical staining of mouse li...



Western blot analysis of extracts of var...