

Product Datasheet

Goat Anti-Mouse IgG (H+L) CY3-conjugated (orb1090778)

Description	Goat Anti-Mouse IgG (H+L) CY3-conjugated
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FACS, ICC, IF, WB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.