

## Product Datasheet

### Goat anti-rat IgG (H+L), HRP-linked Antibody (orb1090781)

<b>Description</b>	Goat anti-rat IgG (H+L), HRP-linked Antibody
<b>Species/Host</b>	Goat
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Note</b>	For research use only
<b>Clonality</b>	Polyclonal
<b>Expiration Date</b>	12 months from date of receipt.