

Product Datasheet

Phospho-NF2 / Merlin (Ser13) Antibody (orb1093974)

Description Phospho-NF2 / Merlin (Ser13) Antibody

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, ICC, IF, IHC-P, WB

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Clonality Polyclonal

Expiration Date 12 months from date of receipt.