
Product Datasheet

Influenza A, H1N1, Hemagglutinin antibody (orb109453)

Description

Mouse monoclonal to Influenza A,

Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Purified influenza virus type A strain H1N1, A/New Caledonia/20/99 strain.
Preservatives	Supplied as a liquid in PBS, pH 7.4, 0.08% sodium azide
Form/Appearance	Supplied as a liquid in PBS, pH 7.4, 0.08% sodium azide
Concentration	See vial for concentration
Storage	May be stored at 4°C for short-term only. For long-term storage and to avoid repeated freezing and thawing, add sterile 40-50% glycerol, aliquot and store at -20°C
Note	For research use only
Application notes	Optimal concentration should be evaluated by serial dilutions.
Isotype	IgG1
Clonality	Monoclonal
Purity	Affinity Purified
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.