


Product Datasheet

Phospho-FXYD1/PLM (Ser83) Antibody (orb1095345)

Description	Phospho-FXYD1/PLM (Ser83) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	A synthesized peptide derived from human FXYD1/PLM around the phosphorylation site of Ser83
Preservatives	Rabbit IgG in phosphate buffered saline, pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol
Concentration	~1mg/ml
Storage	Store at -20 °C. Stable for 12 months from date of receipt
Note	For research use only
Application notes	For western blot, incubate membrane with diluted primary Ab in 5% w/v milk , 1X TBS, 0.1% Tween®20 at 4°C with gentle shaking, overnight
Clonality	Polyclonal
Purity	The Ab is from purified rabbit serum by affinity purification via sequential chromatography on phospho-peptide and non-phospho-peptide affinity columns
Dilution Range	WB 1:500-1:2000, IHC 1:50-1:200
Expiration Date	12 months from date of receipt.



orb1095345
at 1/100
staining Rat
heart t...