

## Product Datasheet

# Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated (orb1095815)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated  |
| <b>Species/Host</b>        | Goat  |
| <b>Reactivity</b>          | Rabbit  |
| <b>Conjugation</b>         | Cy3   |
| <b>Tested Applications</b> | FC, IF  |
| <b>Immunogen</b>           | The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| <b>Preservatives</b>       | 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6   |
| <b>Form/Appearance</b>     | Liquid  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |