

Product Datasheet

Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated (orb1095815)

| | |
|----------------------------|---|
| Description | Goat Anti-Rabbit IgG(H+L) Antibody; Cy3-conjugated |
| Species/Host | Goat |
| Reactivity | Rabbit |
| Conjugation | Cy3 |
| Tested Applications | FC, IF |
| Immunogen | The antibody was purified from antisera by immunoaffinity chromatography using antigens coupled to agarose beads. |
| Preservatives | 0.01M Sodium phosphate, 0.25M NaCl, 50% glycerol, pH7.6 |
| Form/Appearance | Liquid |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)