

Product Datasheet

S Antibody, Biotin conjugated (orb1095818)

Description	S Antibody, Biotin conjugated
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Novel Coronavirus Spike glycoprotein (S) (16-685aa)
Target	S
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	PODTC2
Expiration Date	12 months from date of receipt.