

Product Datasheet

ABHD10 Antibody (orb1095827)

Description	ABHD10 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from human protein . at AA range: 220-300
Target	ABHD10
Preservatives	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9NUJ1
Dilution Range	WB:1:500-2000, ELISA:1:5000-20000
Expiration Date	12 months from date of receipt.