

## Product Datasheet

### Nucleoprotein (NTD) Antibody (orb1095842)

<b>Description</b>	Nucleoprotein (NTD) Antibody
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC
<b>Immunogen</b>	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
<b>Target</b>	N
<b>Preservatives</b>	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P0DTC9</a>
<b>Dilution Range</b>	WB:1000-5000, ELISA:1:1000-1:5000
<b>Expiration Date</b>	12 months from date of receipt.