

Product Datasheet

N Monoclonal Antibody (orb1095843)

Description	N Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa)
Target	N
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	PODTC9
Dilution Range	WB:1000-5000, ELISA:1:1000-1:5000
Expiration Date	12 months from date of receipt.