
Product Datasheet

**Recombinant Human Receptor tyrosine-protein
kinase erbB-2(ERBB2), partial (Active)
(orb1095894)**

Description

Recombinant Human Receptor tyrosine-protein kinase erbB-

Endotoxins	Less than 1.0 EU/ug as determined by LAL method.
Preservatives	Lyophilized from a 0.2 µm filtered 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0
Form/Appearance	Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	C-terminal 6xHis-tagged
Note	For research use only
Application notes	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Protein Sequence	TQVCTGTDMLRLPASPETHLDMLRHLYQGCQVVGQGNLELYLPTNASLSFLQD IQEVQGYVLIAHNQVRQVPLQRLRIVRGTQLFEDNYALAVLDNGDPLNNTTPVTG ASPGGLRELQLRSLTEILKGGVLIQRNPQLCYQDTILWKDIFHKNNQLALTLIDTN RSRACHPCSPMCKGSRGWESSEDCQSLTRTV CAGGCARCKGPLPTDCCHEQ CAAGCTGPKHSDCLACLHFNHSGICELHCPALVTYNTDTFESMPNPEGRYTFGA SCVTACPYNYLSTDVGSCTLVCPLHNQEVTAEDGTQRCEKCSKPCARVCYGLG MEHLREVRAVTSANIQEFAGCKKIFGSLAFLPESFDGDPASNTAPLQPEQLQVFE TLEEITGYLYISAWPDSLPLDSVFNLQVIRGRILHNGAYSLTLQGLGISWLGLRSL RELGSGLALIHNTLHLCFVHTVPWDQLFRNPHQALLHTANRPEDECVGEGLAC HQLCARGHCWGPPTQCVNCSQFLRGQECVEECRVLQGLPREYV NARHCLPC HPECQPQNGSVTCFGPEADQCACAHYKDPPFCVARCPGSKPDL SYMPIWK PDEEGACQPCPINCTHSCVDLDDKGC PAEQRASPLT
Purity	Greater than 95% as determined by SDS-PAGE.
Source	Mammalian cell
MW	70.2 kDa
Uniprot ID	P04626
Expiration Date	6 months from date of receipt.