


Product Datasheet

M13 bacteriophage (g3p) antibody (orb109590)

Description	Antibody specific for the M13 bacteriophage protein coat, amongst others the g3p protein. Provided a...
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, FACS, IHC-P, WB
Immunogen	M13 phage coat proteins.
Concentration	See vial for concentration.
Storage	Product should be stored at -20°C. Aliquot to avoid freeze/thaw cycles
Note	For research use only
Application notes	ELISA Western Blot Flow Cytometry Immunohistochemistry
Isotype	IgG3
Clonality	Monoclonal
Source	Hybridoma produced by the fusion of splenocytes from mice immunized with isolated M13 phage coat proteins and mouse myeloma cells.
Uniprot ID	P69168
Hazard Information	This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals.
Expiration Date	12 months from date of receipt.



Line graph illustrates about the Ag-Ab r...