
Product Datasheet

Recombinant Human P-selectin (SELP), partial (orb1096083)

Description

Recombinant Human P-selectin(SELP),partial

Endotoxins

Not test.

Preservatives

If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.

Form/Appearance

Liquid or Lyophilized powder

Storage

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Tag

C-terminal hFc-Myc-tagged

Note

For research use only

Application notes

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Protein Sequence

WTYHYSTKAYSWNISRKYCQNRYTDLVAIQNKNEIDYLNKVLPPYSSYYWIGIRK
 NNKTWTWVGTKKALTNEAENWADNEPNNKRNNEDCVEIYIKSPSAPGKWNDE
 HCLKKKHALCYTASCQDMSCSKQGELETIGNYTCSCYPGFYGPCEYVRECGE
 LELPQHVLNMCNCSHPLGNFSFNSQCSFHCTDGYQVNGPSKLECLASGIWTKPPQ
 CLAAQCPPLKIPERGNMCLHSAKAFQHQSSCSFSCEEGFALVGPEVVQCTASG
 VWTAPAPVCKAVQCQHLEAPSEGTMDCVHPLTAFAYGSSCKFECQPGYRVRL
 DMLRCIDSGHWSAPLPTCEAISCEPLESPVHGSMDCSPSLRAFQYDTNCSFRCA
 EGFMLRGADIVRCNDNLGQWTAPAPVCQALQCQDLPVPNEARVNCNSHPFGAFRY
 QSVCSFTCNEGLLLVGASVLQCLATGNWNSVPPECQAIPCTPLLSPQNGTMTCV
 QPLGSSSYKSTCQFICDEGYLSGPERLDCTRSGRWTDSPPMCEAIKPELFAPE
 QGSLDCSDTRGEFNVGSTCHFSCDNGFKLEGPNNECTTSGRWSATPPTCKGI
 ASLPTPGLQCPALTPGQGTMYCRHHPGTFGFNTTCYFGCNAGFTLIGDSTLSCR
 PSGQWTAVTPACRAVKCSELHVNKPIAMNCSNLWGNFSYGSICSFHCLEGQLLN
 GSAQTACQENGHWSTTVPTCQAGPLTIQEA

Purity

Greater than 90% as determined by SDS-PAGE.

Source

Mammalian cell

MW

108.8 kDa

Uniprot ID
[P16109](#)
Expiration Date

6 months from date of receipt.