

---

## Product Datasheet

**Recombinant Human Novel Coronavirus  
Nucleoprotein(N)(R203K,G204R) (Active)  
(orb1096702)**

**Description**

Recombinant Human Novel Coronavirus

**Endotoxins**

Not test.

**Preservatives**

Lyophilized from 20 mM Tris-HCl, 0.5 M NaCl, 6% Trehalose, pH 8.0

**Form/Appearance**

Lyophilized powder

**Storage**

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

**Tag**

N-terminal 6xHis-tagged

**Note**

For research use only

**Application notes**

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

**Protein Sequence**

MSDNGPQNQRNAPRITFGGSPDSTGSNQNGERSGARSKQRRPQGLPNNTASW  
FTALTQHGKEDLKFPRGQVPINTNSSPDDQIGYYRRATRRIRGGDGKMKDLSP  
RWYFYLLGTGPEAGLPYGANKDGIWVATEGALNTPKDHIGTRNPANNAIVLQL  
PQGTTLPGGFYAEGSRGGSQASSRSSRSRNSSRNSTPGSSKRTSPARMAGNG  
GDAALALLLLDRLNQLESKMSGKQQQQGQTVTKKSAEASKKPRQKRTATKA  
YNVTQAFGRRGPEQTQGNFGDQELIRQGTDYKHWPQIAQFAPSASAFFGMSRI  
GMEVTPSGTWLTYTGAIKLDDKDPNFKDQVILLNKHIDAYKTFPTEPKKDKKKK  
ADETQALPQRQKKQQTVLLPAADLDDFSKQLQSMSSADSTQA

**Purity**

Greater than 96% as determined by SDS-PAGE.

**Source**

E.coli

**MW**

49.8 kDa

**Uniprot ID**[P0DTC9](#)**Expiration Date**

6 months from date of receipt.