

Product Datasheet

Connexin 43 (phospho Ser261) Polyclonal Antibody (orb1097299)

Description	Connexin 43 (phospho Ser261) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P17302
Expiration Date	12 months from date of receipt.