

Product Datasheet

Ephrin-B1/2/3 (phospho Tyr324) Polyclonal Antibody (orb1097364)

Description	Ephrin-B1/2/3 (phospho Tyr324) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q15768, P52799, P98172
Expiration Date	12 months from date of receipt.