

Product Datasheet

KCNQ2/3/4/5 (phospho Thr217/246/223/251) Polyclonal Antibody (orb1097489)

Description	KCNQ2/3/4/5 (phospho Thr217/246/223/251) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O43526 , Q9NR82 , O43525 , P56696
Expiration Date	12 months from date of receipt.