

Product Datasheet

Ku-80 (phospho Thr714) Polyclonal Antibody (orb1097496)

Description	Ku-80 (phospho Thr714) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P13010
Expiration Date	12 months from date of receipt.