

Product Datasheet

Lyn (phospho-Tyr507) rabbit pAb Antibody (orb1097516)

Description	Lyn (phospho-Tyr507) rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P07948
Expiration Date	12 months from date of receipt.