

Product Datasheet

MAPKAPK-2 (phospho Thr222) Polyclonal Antibody (orb1097520)

Description	MAPKAPK-2 (phospho Thr222) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P49137
Expiration Date	12 months from date of receipt.