

Product Datasheet

MARK1/2/3/4 (phospho Thr215) Polyclonal Antibody (orb1097525)

| | |
|----------------------------|---|
| Description | MARK1/2/3/4 (phospho Thr215) Polyclonal Antibody |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | ELISA, IHC-P |
| Storage | Store at -20°C for long term preservation (recommended). Store at 4°C for short term use. |
| Note | For research use only |
| Clonality | Recombinant |
| Uniprot ID | P27448 , Q9P0L2 , Q96L34 , Q7KZ17 |
| Expiration Date | 12 months from date of receipt. |