

Product Datasheet

MKK7 (phospho-Ser271/Thr275) rabbit pAb Antibody (orb1097549)

Description	MKK7 (phospho-Ser271/Thr275) rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O14733
Expiration Date	12 months from date of receipt.