

Product Datasheet

Moesin/Ezrin/Radixin (phospho Thr558) Polyclonal Antibody (orb1097555)

Description	Moesin/Ezrin/Radixin (phospho Thr558) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P15311 , P35241 , P26038
Expiration Date	12 months from date of receipt.