

Product Datasheet

NFATc4 (phospho Ser676) Polyclonal Antibody (orb1097582)

Description	NFATc4 (phospho Ser676) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q14934
Expiration Date	12 months from date of receipt.