

## Product Datasheet

# PDE4B/C/D (phospho Ser133/119/190) Polyclonal Antibody (orb1097630)

<b>Description</b>	PDE4B/C/D (phospho Ser133/119/190) Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<a href="#">Q07343</a> , <a href="#">Q08493</a> , <a href="#">Q08499</a>
<b>Expiration Date</b>	12 months from date of receipt.