

Product Datasheet

PKC (phospho Thr497) Polyclonal Antibody (orb1097638)

Description	PKC (phospho Thr497) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P05771 , Q04759 , Q05513 , P17252 , P24723 , Q02156 , P05129 , Q05655
Expiration Date	12 months from date of receipt.