

Product Datasheet

PNK1 (phospho-Ser114/Thr118) rabbit pAb Antibody (orb1097652)

Description	PNK1 (phospho-Ser114/Thr118) rabbit pAb Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q96T60
Expiration Date	12 months from date of receipt.