

Product Datasheet

RGS16 (phospho Tyr168) Polyclonal Antibody (orb1097692)

Description	RGS16 (phospho Tyr168) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	O15492
Expiration Date	12 months from date of receipt.