

## Product Datasheet

# Ribosomal Protein S6 (phospho Ser235/S236) Polyclonal Antibody (orb1097696)

<b>Description</b>	Ribosomal Protein S6 (phospho Ser235/S236) Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IHC-P, WB
<b>Storage</b>	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
<b>Note</b>	For research use only
<b>Clonality</b>	Recombinant
<b>Uniprot ID</b>	<b>P62753</b>
<b>Expiration Date</b>	12 months from date of receipt.