

Product Datasheet

Ribosomal Protein S6 (phospho Ser240) Polyclonal Antibody (orb1097697)

Description	Ribosomal Protein S6 (phospho Ser240) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	P62753
Expiration Date	12 months from date of receipt.