

Product Datasheet

Rrn3 (phospho Ser649) Polyclonal Antibody (orb1097712)

Description	Rrn3 (phospho Ser649) Polyclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.
Note	For research use only
Clonality	Recombinant
Uniprot ID	Q9NYV6
Expiration Date	12 months from date of receipt.