

## Product Datasheet

### Exendin 4 (3-39) peptide (orb1097933)

<b>Description</b>	Exendin 4 (3-39) peptide
<b>Form/Appearance</b>	White lyophilized powder
<b>Storage</b>	Keep in cool and dry place
<b>Note</b>	For research use only
<b>Protein Sequence</b>	Glu-Gly-Thr-Phe-Thr-Ser-Asp-Leu-Ser-Lys-Gln-Met-Glu-Glu-Glu-Ala-Val-Arg-Leu-Phe-Ile-Glu-Trp-Leu-Lys-Asn-Gly-Gly-Pro-Ser-Ser-Gly-Ala-Pro-Pro-Pro-Ser-NH <sub>2</sub>
<b>Purity</b>	> 95%
<b>Expiration Date</b>	12 months from date of receipt.