



Product Datasheet

Claudin18/CLDN18 Antibody (orb1098020)


Description	Claudin18/CLDN18 Antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human Claudin18/CLDN18 recombinant protein (Position: S31-V261).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	P56856
Expiration Date	12 months from date of receipt.



IHC analysis of Claudin18/CLDN18 using a...



IHC analysis of Claudin18/CLDN18 using a...



IF analysis of Claudin18/CLDN18 using an...