



Product Datasheet

Tryptophan Hydroxylase/TPH1 Antibody (orb1098051)


Description	Tryptophan Hydroxylase/TPH1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IF, IHC, WB
Immunogen	E.coli-derived human Tryptophan Hydroxylase/TPH1 recombinant protein (Position: K383-I444).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	60 kDa
Uniprot ID	P17752
Expiration Date	12 months from date of receipt.



IHC analysis of Tryptophan Hydroxylase/T...



IHC analysis of Tryptophan Hydroxylase/T...



IHC analysis of Tryptophan Hydroxylase/T...