



Product Datasheet

CSNK2A2 Antibody (orb1098102)


Description	CSNK2A2 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, WB
Immunogen	A synthetic peptide corresponding to a sequence of human CSNK2A2 (EHPYFYPVVKEQSQPCADNAVLSSGLTAAR).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	P19784
Expiration Date	12 months from date of receipt.



IHC analysis of CSNK2K2 using anti-CSNK2...



IHC analysis of CSNK2K2 using anti-CSNK2...



IHC analysis of CSNK2K2 using anti-CSNK2...