


Product Datasheet

Human CCL14 ELISA Kit (DIY Antibody Pairs) (orb1098126)

Description	Human CCL14 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: T20-N93
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA)
Uniprot ID	Q16627
Expiration Date	Please enquire.



Human
CCL14 EZ
Set ELISA
Kit
standard
CU...