


Product Datasheet

Human Decorin / DCN ELISA Kit (DIY Antibody Pairs) (orb1098129)

Description	Human Decorin / DCN ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D31-K359
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39747 MW
Uniprot ID	P07585
Expiration Date	Please enquire.



Human Decorin EZ Set ELISA Kit standard ...