


Product Datasheet

Human IL-8 ELISA Kit (DIY Antibody Pairs) (orb1098134)

Description	Human IL-8 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml (cell culture supernates) 7.8 pg/ml - 500 pg/ml (serum, plasma)
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S28-S99
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum and plasma (heparin, EDTA, citrate).
MW	11098 MW
Uniprot ID	P10145
Expiration Date	Please enquire.



Human IL-8
EZ Set
ELISA Kit
standard
cur...