


Product Datasheet

Human CRP ELISA Kit (DIY Antibody Pairs) (orb1098150)

Description	Human CRP ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: F17-P224
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	P02741
Expiration Date	Please enquire.



Human CRP
EZ-Set
ELISA Kit
standard
curv...