


Product Datasheet

Human Cystatin C ELISA Kit (DIY Antibody Pairs) (orb1098151)

Description	Human Cystatin C ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: M1-A146
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, urine and human milk.
Uniprot ID	P01034
Expiration Date	Please enquire.



Human
Cystatin C
EZ-Set
ELISA Kit
standa...