


Product Datasheet

Human IL-33 ELISA Kit (DIY Antibody Pairs) (orb1098153)

Description	Human IL-33 ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: S112-T270
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) and cell lysates.
Uniprot ID	095760
Expiration Date	Please enquire.



Human IL-33 EZ-Set ELISA Kit standard CU...