


Product Datasheet

Human APOE ELISA Kit (DIY Antibody Pairs) (orb1098157)

Description	Human APOE ELISA Kit (DIY Antibody Pairs)
Reactivity	Human
Range	3.12 ng/ml - 200 ng/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: K19-H317
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva, milk and urine.
Uniprot ID	P02649
Expiration Date	Please enquire.



Human APOE EZ-Set ELISA Kit standard cur...