


Product Datasheet

Mouse CCL20/MIP-3 alpha ELISA Kit (DIY Antibody Pairs) (orb1098162)

Description	Mouse CCL20/MIP-3 alpha ELISA Kit (DIY Antibody Pairs)
Reactivity	Mouse
Range	7.8 pg/ml - 500 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A28-M97
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
Uniprot ID	089093
Expiration Date	Please enquire.



Mouse CCL20/MIP-3 alpha EZ-Set ELISA Kit...